

746 Data File Format Definition

Subject: 746 Data Format
Date: 7 April 1999
Author: John Cooper
Status: Final

The 746 program, 746SIG-B.EXE, generates two data files: .ALL and .AVG. Generation of the .AVG file is an operator selectable option. These files contain comments and data pertaining to a source scan.

Files are ASCII, space delimited, CR/LF terminated. The program runs under DOS, thus the 8.3 filename limitation applies.

Filenames: SYMMDDA.ext
 where S

The source being scanned:

A - Ames 30 inch Sphere
B - Boreas Hemisphere
C - CHORS Sphere
F - Flat Plate (any type)
H - GSFC 48-in Hemisphere
J - MISR Sphere at JPL or ASTER 1M Sphere at NEC
K - MELCO 1M SWIR Sphere
L - Standard Lamp (Halogen type, any model)
M - Ames 20 inch Hemisphere
R - SIS 18-in aperture Sphere
S - GSFC 6-ft. Sphere
T - GSFC 12-in Sphere
X - Test scan

YY Two digit year
MM Two digit month
DD Two digit day
A An alphanumeric showing the chronological position of the file
for that day, A indicates first scan, B indicates second, ...
ext File contents indicator:
ALL Raw data file, containing all sampled data points
for each wavelength
AVG Averaged data file, only the averaged raw data for
each wavelength

Raw data (ALL) Format

There are two sections in this type of file: a five line header and the data.

Italics indicate variable information.

ALL Header, in order:

<i>desg/com</i>	Source designation and comments, variable length
System output [<i>units</i>]	Definition of data units. <i>units</i> can be Volts or Amps
Date: YYMMDD	Date of data collection
<i>snm</i> to <i>enm</i> <i>inm</i> steps <i>n</i> samples	<i>s</i> - Starting wavelength of scan, in nanometers <i>e</i> - Ending wavelength of scan, in nanometers <i>i</i> - Interval between sampled wavelengths, in nanometers <i>n</i> - Number of samples taken if noise limits are not exceeded. If the noise limits are exceeded, more samples are taken, up to the system maximum of 100 samples.
<i>tC rh%RH distcm</i> dist <i>smm</i> slits	<i>t</i> - Temperature entered by technician, in °C <i>rh</i> - Relative Humidity entered by technician, in % <i>dist</i> - Distance between instrument and source, in cm <i>s</i> - Spectroradiometer slit size, in mm

ALL data:

There are as many data lines as there are sampled wavelengths. The data line has the following structure:

λ hh:mm:ss s_1 s_2 s_3 ... s_n s_{n+1} s_{n+2} s_m

where: λ sampled wavelength

hh:mm:ss	time at start of sampling
s ₁	first sample
s ₂	second sample
s _n	n th sample
s _{n+1}	n+1 th sample, taken only if noise limits are exceeded
s _{n+2}	n+2 th sample, taken only if noise limits are exceeded
s _m	m th sample, taken only if noise limits are exceeded. m _{max} is 100.

Sample ALL file:

Filename is L990326b.all

```
f463 blocked
System output [Volts]
Date: 990326
800nm to 2400nm 20nm steps 5 samples
24c 30%RH 50cm dist. 2.5nm slits
800 17:15:40 5.000E-06 4.300E-06 3.000E-06 0.000E+00 1.200E-06 ... 4.500E-06
820 17:15:52 2.800E-06 4.000E-06 5.200E-06 1.700E-06 4.000E-06 ... 9.000E-07
...
2380 17:32:08 2.900E-06 5.500E-06 3.300E-06 2.200E-06 2.900E-06 ... 1.000E-05
2400 17:32:19 1.200E-06 6.700E-06 5.000E-07 9.000E-07 4.400E-06 ... 4.000E-06
```

Averaged data (AVG) format

There are two sections in this type of file: a one line header and the data.

Italics indicate variable information.

AVG Header:

The averaged file header is a single line in the following format:

"*desg/com*", "System output [*units*]", *YYMMDD,s,e,i*

where

<i>desg/com</i>	Source designation and comments, variable length
System output [<i>units</i>]	Definition of data units. <i>units</i> can be Volts or Amps
<i>YYMMDD</i>	Date of data collection
<i>s</i>	Starting wavelength of scan, in nanometers
<i>e</i>	Ending wavelength of scan, in nanometers
<i>i</i>	Interval between sampled wavelengths, in nanometers

AVG data:

There are as many data lines as there are sampled wavelengths. A data line is a single number representing the averaged data.

Sample AVG file:

Filename is L990326b.avg

```
"f463 blocked", "System output [Volts]", 990326, 800, 2400, 20
2.365000E-06
2.615000E-06
3.445000E-06
3.635000E-06
...
1.143500E-05
1.042500E-05
8.500000E-06
7.855000E-06
```